

Labnet Mini Microcentrifuges

- Quick release rotor system
- Electronic braking
- Dual lid switch
- Compact, rugged design
- Each centrifuge box includes two rotors, three tube adapters and a power adapter



SPECIFICATIONS

| | |
|------------------|--|
| Maximum Speed | 6,000 rpm / 2,000 x g |
| Maximum Capacity | 8 x 1.5/2.0 mL or 4 x 0.2mL strips |
| Dimensions | 6.8 x 5.8 x 4.8 in / 17.2 x 14.7 x 12.2 cm |
| Electrical | AC 100 - 240V, 50/60Hz |

ORDERING INFORMATION

| | |
|-------------|--|
| C1601-B | Labnet Mini Microcentrifuge, blue lid, 100-240V |
| C1601-P | Labnet Mini Microcentrifuge, purple lid, 100-240V |
| C1601-R | Labnet Mini Microcentrifuge, red lid, 100-240V |
| C1601-G | Labnet Mini Microcentrifuge, gray lid, 100-240V |
| 6770-RT | Rotor for 8 x 1.5/2.0mL tubes |
| 6770-RTS | Rotor for 4 x 8-strip tubes |
| 6770-RTA.20 | Individual adapters for 0.2mL tubes, fits in rotor 6770-RT, pack of 8 |
| 6770-RTA.25 | Individual adapters for 0.25mL tubes, fits in rotor 6770-RT, pack of 8 |
| 6770-RTA.5 | Individual adapters for 0.5mL tubes, fits in rotor 6770-RT, pack of 8 |
| 401025 | Power adapter, 24V |