

LABNET INTERNATIONAL, INC.

**BIOSALAB**  
MOLECULAR BIOLOGY SOLUTIONS

# MPS 1000™ Mini Plate Spinner

Finally, a compact and economical centrifuge for PCR plates

**Labnet**  
Labnet International, Inc.



- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges

Patent Pending



## SPECIFICATIONS

Capacity:	2 PCR Plates
Rotor Speed:	2500 rpm
G-Force:	500 xg
Rotor:	Fixed Vertically
Dimensions (WxDxH):	7.5 x 8.2 x 7.2 in. 19 x 21 x 18.3 cm
Weight:	6 lbs / 3 kg
Electrical:	120V or 230V, 50-60 Hz

## ORDERING INFORMATION

Catalog No.	Description
C1000*	<b>Mini Plate Spinner</b> centrifuge with two position rotor

\* To order this product in 230V, please add (-230V) to the item number.

**Labnet**  
Labnet International, Inc.

PO Box 841 • Woodbridge • NJ • 07095 • USA • website: <http://www.labnetlink.com>  
US • toll free: 888-LABNET1 • fax: (732) 417-1750 • email: [labnet@labnetlink.com](mailto:labnet@labnetlink.com)  
International: +1-732-417-0700 • fax: +1-732-417-1750 • [international@labnetlink.com](mailto:international@labnetlink.com)